

3M Science.
Applied to Life.™

3M™ Impregum™ Super Quick
Polyether Impression Material

**Capture every
detail in 2 minutes?
Yes, it can.**



**Now also
available in
cartridge
delivery!**

The fastest way to polyether smaller cases. 3M™ Impregu

Impregum Polyether Impression Material has always been the trusted solution for complex situations and implants. To ensure that the highest demands are fulfilled without limiting productivity, 3M Oral Care created Impregum Super Quick Polyether Impression Material.

It combines the outstanding performance and reliability of a polyether with the speed of a VPS. Deep margins and moist environments will no longer slow you down.



Quicker set.

So fast, so good: With just a two minute setting time Impregum Super Quick Material allows for faster impressions in less complex cases like single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units.



100% Polyether. Trusted precision.

Like all polyethers Impregum Super Quick Materials demonstrate excellent performance in moist situations and great flow behavior. Both lead to a high level of accuracy and detail reproduction for high-precision restorations.

It is especially important that you get the full working time of your impression material. Because it is a polyether, you get the full 45 seconds of working time with Impregum Super Quick Material. You are able to avoid inaccuracies associated with premature setting that is experienced with VPS materials.



Make your patients feel comfortable.

An improved fresh, minty taste together with a short time in the mouth significantly increases your patients' comfort.



Polyether meets CAD/CAM workflow.

Impregum Super Quick Material is scannable with the current generation of scanners for a more efficient cooperation with your lab.

precision for m™ Super Quick.



Fast and gentle moisture management and retraction of gingival tissue with 3M™ Astringent Retraction Paste.




Convenient and highly accurate application of 3M™ Impregum™ Super Quick Light Body Polyether Impression Material using 3M™ Intra-oral Syringe Purple.



Precise, void-free details plus contrasting colors for excellent readability.

High-precision 1-step impression made with 3M™ Impregum™ Super Quick Polyether Impression Material (Heavy and Light Body) and single-use 3M™ Dual Arch Impression Tray.



3M™ Impregum™ Super Quick


HIGHLY SATISFIED






A clinical evaluation shows a **very high satisfaction rate** of almost **95%** regarding the **precision** of 3M™ Impregum™ Super Quick Material impressions and the **fit of the final restorations**.



Source: 3M field evaluation with 447 participants from Europe and the U.S., 2017.

- First superfast polyether with 45 sec. working time and quick 2-minute setting
- Available in heavy, medium, light and ultra light body consistencies
- For monophase and 1-step impression techniques
- Perfect for impressions for single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units



	Dispensing system	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
3M™ Impregum™ Penta™ Super Quick (Medium Body)		Monophase or 1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Medium Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Ultra Light Body)		1-step	✓	✓	✓	✓	0:45	2:00

Ordering information



3M™ Impregum™ Super Quick Medium Body/Light Body/Ultra Light Body Introductory Kit (69421)



3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body Introductory Kit (69415)



3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body Introductory Kit (69407)

Item no. Product information

Introductory Kit for 3M™ Garant™ Dispenser

69421 3M™ Impregum™ Super Quick Medium Body/Light Body/Ultra Light Body
2 3M™ Impregum™ Super Quick MB cartridges – 50 ml each, 1 3M™ Impregum™ Super Quick LB cartridge – 50 ml, 1 3M™ Impregum™ Super Quick ULB cartridge – 50 ml, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive

Introductory Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit

69414 3M™ Impregum™ Penta™ Super Quick Medium Body
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 1 3M™ Penta™ Elastomer syringe, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

69415 3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

69407 3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body
1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ Super Quick HB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

Refills

69385 3M™ Impregum™ Penta™ Super Quick Medium Body
2 base pastes – 300 ml each, 2 catalysts – 60 ml each

69417 3M™ Impregum™ Super Quick Medium Body
4 cartridges – 50 ml each, 10 3M™ Garant™ Mixing Tips purple

69406 3M™ Impregum™ Penta™ Super Quick Heavy Body
2 base pastes – 300 ml each, 2 catalysts – 60 ml each

69379 3M™ Impregum™ Super Quick Light Body
4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple

69418 3M™ Impregum™ Super Quick Ultra Light Body
4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple

Accessories

69387 3M™ Impregum™ Penta™ Super Quick Medium Body Cartridge

69401 3M™ Impregum™ Penta™ Super Quick Heavy Body Cartridge

71622 48 3M™ Dual Arch Impression Trays Posterior

71623 35 3M™ Dual Arch Impression Trays Anterior

30600 3M™ Polyether Tray Adhesive
1 bottle – 17 ml



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-634-2249

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

3M.com/Dental

3M, Garant, Impregum, Penta and Pentamix are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies.

© 3M 2020. All rights reserved. 70-2013-7072-6