

3M Science.
Applied to Life.™

3M™ Impregum™ Super Quick Polyether Impression Material

Capture every
detail in 2 minutes?
Yes, it can.




Now also
available in
cartridge
delivery!



Quicker set.

With just a two minute setting time, Impregum Super Quick Material lives up to its name and allows for faster impressions in cases like single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units.



100 % Polyether. Trusted precision.

Impregum Super Quick Materials demonstrate excellent performance in moist situations and exhibit great flow behavior, leading to higher levels of accuracy and detailed, high-precision restorations.

The nature of polyether allows you to get the full 45 seconds of working time, avoiding inaccuracies associated with premature setting times common with VPS materials.




Make your patients feel comfortable.

An improved fresh, minty taste together with a short time in the mouth significantly increases your patients' comfort.

3M™ Impregum™ Super Quick Polyether Impression Material

- For monophasic and 1-step impression techniques
- Perfect for impressions for single-unit crowns, implants, inlays/onlays and smaller bridges up to 3 units
- Available in heavy, medium, light and ultra light body consistencies
- Impregum Super Quick Material is scannable with the current generation of scanners for a more efficient workflow.



The fastest way to polyether precision for smaller cases.

	Dispensing system	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
3M™ Impregum™ Penta™ Super Quick (Medium Body)		Monophase or 1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Medium Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)		1-step	✓	✓	✓	✓	0:45	2:00
3M™ Impregum™ Super Quick (Ultra Light Body)		1-step	✓	✓	✓	✓	0:45	2:00

Ordering information



3M™ Impregum™ Super Quick Medium Body/Light Body/Ultra Light Body Introductory Kit (69421)



3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body Introductory Kit (69415)



3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body Introductory Kit (69407)

Item no. Product information

Introductory Kit for 3M™ Garant™ Dispenser

69421 3M™ Impregum™ Super Quick Medium Body/Light Body/Ultra Light Body
2 3M™ Impregum™ Super Quick MB cartridges – 50 ml each, 1 3M™ Impregum™ Super Quick LB cartridge – 50 ml, 1 3M™ Impregum™ Super Quick ULB cartridge – 50 ml, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive

Introductory Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit

69414 3M™ Impregum™ Penta™ Super Quick Medium Body
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 1 3M™ Penta™ Elastomer syringe, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

69415 3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ Super Quick MB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

69407 3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body
1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ Super Quick HB cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes – Purple, 5 3M™ Garant™ Mixing Tips – Purple, 5 3M™ Garant™ Intra-oral Tips – White, 1 3M™ Polyether Tray Adhesive, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Dual Arch Trays (3 posterior, 3 anterior trays)

Refills

69385 3M™ Impregum™ Penta™ Super Quick Medium Body
2 base pastes – 300 ml each, 2 catalysts – 60 ml each

69417 3M™ Impregum™ Super Quick Medium Body
4 cartridges – 50 ml each, 10 3M™ Garant™ Mixing Tips purple

69406 3M™ Impregum™ Penta™ Super Quick Heavy Body
2 base pastes – 300 ml each, 2 catalysts – 60 ml each

69379 3M™ Impregum™ Super Quick Light Body
4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple

69418 3M™ Impregum™ Super Quick Ultra Light Body
4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple

Accessories

69387 3M™ Impregum™ Penta™ Super Quick Medium Body Cartridge

69401 3M™ Impregum™ Penta™ Super Quick Heavy Body Cartridge

71622 48 3M™ Dual Arch Impression Trays Posterior

71623 35 3M™ Dual Arch Impression Trays Anterior

30600 3M™ Polyether Tray Adhesive
1 bottle – 17 ml



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-634-2249

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

[3M.com/Dental](https://www.3m.com/Dental)

3M, Garant, Impregum, Penta and Pentamix are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies.
© 3M 2020. All rights reserved. 70-2013-7071-8