



# SafeMask<sup>®</sup> Architect Pro<sup>™</sup>

Surgical N95 Respirator



PROUDLY MADE IN CANADA  
**BY THE MEDICAL  
MASK EXPERTS**


# Reliable respiratory protection that filters out over 95% of particulate aerosols free of oil

Adjustable nosepiece

Unique single-split comfort bands

Lightweight, breathable construction

High-quality materials and duck-bill shape



Made according to recognized biocompatibility, flammability and fluid resistance standards for use in healthcare facilities, the Medicom SafeMask® Architect Pro™ ensures reliable respiratory protection by filtering more than 95% of airborne particles free of oil. Its innovative “duck bill” shape provides optimal breathability with no exhalation valve required. Provides up to 8-hours of protection. Available in 3 sizes.

Medicom SafeMask Architect Pro N95 surgical respirator has full approval from the National Institute for Occupational Safety and Health (NIOSH) of the Centers for Disease Control and Prevention (CDC) under approval number TC-84A-9375. NIOSH certification means that the respirator has met the requirements of Title 42, Code of Federal Regulations, Part 84 (42 CFR 84).



#### Item codes:

**203214** Small      **203314** Medium      **203414** Large  
50 masks/box; 10 boxes/case



Like you, we never compromise. We use carefully selected materials and engineer our products to offer protection you can count on. We have been providing peace of mind to healthcare professionals for over 30 years. We take pride in the protection we offer, so you can take pride in yours.

For more information, please visit [medicom.com](https://www.medicom.com) or contact us at: [cservice@medicom.ca](mailto:cservice@medicom.ca) / 1-800-361-2862

AMD Medicom Inc. 2555 Chemin de l'Aviation, Pointe-Claire, Montreal, Quebec H9P 2Z2 Canada