



Ivory[®] **ReLeaf**[™]

Dentistry's best hands-free suction solution.

An innovative, hands-free HVE suction device that conveniently connects to existing dental vacuum systems. ReLeaf assists with evacuation, retraction and maintaining a dry field during various dental and hygiene procedures.

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

Ivory® ReLeaf™

ReLeaf ensures that every practitioner can continually provide comfort and efficiency to any dental appointment. With 280 degrees of suction, you never have to worry about suction location again.

- Hands-free, high-volume efficiency
- Soft, comfortable and universal mouthpiece for your patients
- Simple setup—connects to existing HVE system in < 60 seconds



Patient: Hands-free comfort and control

- Latex free and BPA free, custom blend polymer enables Leaf to be soft and comfortable
- Patient easily controls fluid buildup and removal, minimizing risk of choking and gagging
- Allows for effortless movement and communication without affecting suction

Practitioner: Ergonomic efficiency for a wide range of procedures

- Maintain effective suction and dry field during wide array of dental and hygiene procedures
- Work in all 4 quadrants with ergonomic comfort and full view of oral cavity - without need to maneuver
- Improves clinical outcomes with greater focus on the details that matter

Practice: Time savings and practice productivity

- Up to 15 minutes of time savings per patient, translating into hundreds of dollars in production throughout the day
- Frees staff capacity allowing more time to be devoted to other aspects of dental treatment
- Significantly more cost-effective than most conventional suction solutions

Patient and Staff Safety

ReLeaf allows for safer dental and hygiene procedures by minimizing aerosol risk and risk of backflow

Aerosol production raises the risk of airborne contaminants. As recommended by the ADA and CDC, the use of high volume evacuation has been shown as the best way to reduce aerosols.

Risk is also involved when patients close their lips around a saliva ejector, as backflow from low-volume saliva ejectors can occur when the pressure in the patient's mouth is less than that in the evacuator. This backflow can include previous patients' waste material, colonized biofilm from the tubing, and chemicals from the tubing, and chemicals from flushing the line is not an issue with ReLeaf.



View From A Patient

"I found ReLeaf to be very comfortable and liked that I was able to communicate with ReLeaf in my mouth. I also liked that I had control over the water and fluid buildup as I was able to close and remove all fluid whenever I wanted to."

For all patients (8 years and above) including:

- TMJ patients
- Special needs patients

View From A Hygienist

"Whether it be prophylaxis, ultrasonic scaling or debridement, ReLeaf lets me relax and focus on the job at hand. Before ReLeaf, I would spend as much time maneuvering the suction as I did actually scaling teeth. Now, ReLeaf basically relieves me of having to worry about and deal with suction. Any product that can save me time and help me provide exceptional care to my patients will get my attention and respect – and that's exactly what ReLeaf does."

Gwen Smukowski, RDH

Hygiene applications & procedures:

- Ultrasonic Scaling
- Root Planing
- Prophylaxis
- Debridements
- Sealants
- Whitening


View From A Dentist

"I appreciate the fact that ReLeaf can save up to 15 minutes of time per procedure. When you multiply that by the number of hygiene procedures we do each day, that has a huge impact on practice productivity. And when you factor in the impact of the improved patient experiences on patient retention and recruitment, ReLeaf can really be a wonderful asset for the practice."

Lee Ann Brady, DMD


Dental applications & procedures:

- Restorative Procedures
- Endodontics and Periodontics
- Intraoral Scanning
- Placing Ortho/Aligner Attachments
- Bite Blocks, Rubber Dams, and Dry Angles



Simple Quick Setup

- Less than 60 seconds to set up
- Conveniently connects to existing dental HVE systems and adapts to all standard dental unit hardware
- No extensive training necessary



Ivory® ReLeaf™	
66075335	ReLeaf Standard Kit 100 Leaves (2 x 50 Packs), 2 Long HVE Hoses, 6 Short Quick Disconnect Hoses, 6 U-Shaped Connectors
66075336	ReLeaf Pro Kit 200 Leaves (4 x 50 Packs), 2 Long HVE Hoses, 6 Short Quick Disconnect Hoses, 6 U-Shaped Connectors
66075339	ReLeaf 50 Refill 50 Leaves
66075340	ReLeaf 100 Refill 100 Leaves
66075344	ReLeaf Hose Accessory Kit 2 Long HVE Hoses, 4 Short Quick Disconnect Hoses, 4 U-Shaped Connectors
66075376	ReLeaf U-Shaped Connectors 6 U-Shaped Connector Refills

Distributed by:
Kulzer, LLC
 4315 S. Lafayette Blvd.
 South Bend, IN 46614
 (800) 431-1785