

Purelife

Earloop Procedure Masks, Non-Sterile

Masks meet current ASTM 2100 performance criteria for low barrier mask.

Item Codes	1150101, 1150102
Country of Origin	China
Ear loop Length	165 mm
Nose piece length	114 mm
Material	polypropylene (inner, filter, and outer layer)
Bacterial Filtration Efficiency (BFE)	≥ 95%
Differential Pressure	less than 4
Fluid Resistance	Resistant at 80 mmHg
Textile Flammability	No flame spread before 3.5 sec (Class 1)
Cytotoxicity	Grade = 0
Skin Irritation	Negligible
Sensitization	Negligible
Colors	Blue and Pink
Body Size	Length: 178 mm, Width: 89 mm
Usage	For Single-use only
Packaging	50 per box / 10 boxes per case / 500 masks per case