

TECHNICAL DATA SHEET

Non-Woven Sponges

Intended Use

A hydrophilic non-sterile sponge intended to cover a wound, such as to be placed directly on a patient's wound to absorb exudates.

Characteristic	Specification		
Material Composition	30% Polyester and 70% Viscose		
Description	Soft non-woven material, folded (Refer to Appendix 1 specific attribute)		
Absorbency	Retention ≥ 10 g/g		
Color	White		
Ply	4 ply		
Dimension	Size	Width (cm)	Length (cm)
	2" x 2" (5 x 5 cm)	$\geq 4,9$	$\geq 4,9$
	4" x 4" (10 x 10 cm)	$\geq 9,9$	$\geq 9,9$

Additional Information

Storage Condition	Country of Origin
Store under cool dry place	China

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

PureLife, LLC

Item Code	Description	Packaging Information	Size	Treatment
Outer Case: 1110301C Inner Box: 1110301	PureLife Non-Woven Sponges	200 Sponges/Sleeve, 25 Sleeves /Case, 5000 Sponges/Case Sleeve made of Kraft paper Shipper case is made of corrugated fiberboard	2" x 2" (5 x 5 cm)	Non-Sterile
Outer Case: 1110302C Inner Box: 1110302	PureLife Non-Woven Sponges	200 Sponges/Sleeve, 10 Sleeves /Case, 2000 Sponges/Case Sleeve made of Kraft paper Shipper case is made of corrugated fiberboard	4" x 4" (10 x 10 cm)	Non-Sterile