



SafeMask® Classics™

Procedure Earloop Face Masks



Meet all your daily needs with a single line of masks.

You need reliable respiratory protection that is breathable and comfortable, even during the longest procedures. That's why every SafeMask® Classics™ earloop face mask is made of quality materials that provide the ASTM level protection you need without the gaping, discomfort or irritation often caused by lower-quality masks.

Available colours:    



VERSATILITY TO MEET ALL YOUR DAILY NEEDS



Shingle pleat design

Prevents fluid pooling and cross-contamination



High-quality materials

Provide reliable protection without compromising breathability



Adjustable nosepiece

Ensures a comfortable, customized fit



Soft, stretchy earloops

Prevent pulling and irritation

- **Non-woven spunbond polypropylene outer layer** provides maximum protection against fluid penetration.
- **High quality media** ensures optimal filtration and breathability.
- Non-woven spunbond polypropylene **inner layer** is made of dye-free material that will not lint or tear during use.
- Ample length **earloops** are soft, strong, ultrasonically sealed and attached to the outside of the mask to prevent pulling and irritation.
- Adjustable, **enclosed plastic-covered iron nose piece** provides customized protection and comfort.
- **Shingle pleat design** prevents fluid pooling and cross contamination.
- Available in all **3 ASTM levels**.



ASTM
LEVEL 1

205114	White
205115	Blue
205116	Pink
205117	Lavender

ASTM
LEVEL 2

205214	White
205215	Blue
205216	Pink
205217	Lavender

ASTM
LEVEL 3

205314	White
205315	Blue
205316	Pink
205317	Lavender



Like you, we never compromise. We use carefully selected materials and engineer our products to offer protection you can count on. We have been providing peace of mind to healthcare professionals for over 30 years. We take pride in the protection we offer, so you can take pride in yours.

For more information, please visit medicom.com or contact us at: cservice@medicom.ca / 1-800-361-2862

AMD Medicom Inc. 2555 Chemin de l'Aviation, Pointe-Claire, Montreal, Quebec H9P 2Z2 Canada